



FFP2 FACE MASKS

€0.71

SKU: 88360

Categories: [Bioide](#), [FaceProtect](#), [Masks](#),
[Offers](#), [PersonProtect](#)

PRODUCT DESCRIPTION

Disposable FFP2 Face Mask

Material: 4ply(non-woven+Hot air cotton+meltblown cloth+non-woven)

Color: White

Standard: China GB2626-2006

Application:Personal protection,PM2.5 protection

Filtrating effect:95%

Product Features

- 1.Multi layer filter structure protection
- 2.Ear belt / head-mounted
- 3.Comfortable and breathable
- 4.Good protective effect
- 5.Comfortable soft elastic band
- 6.Soft sponge nose pad
- 7.Aluminum adjustable nose clip

General information - Contact us for more details if needed:

Supplied from multiple factories

All masks quality checked and CE certificates verified
CE certificate / Test Certificates / Photos can vary by manufacturer.
Packaging: Variable. Op-Bags 1/4 masks / Larger bags 10/20/50 masks.
Typical QTY in carton: 1000 masks.
Special packing configuration available upon request.